

Product datasheet for **HP230314**

TCTN1 Human qPCR Primer Pair (NM_001082538)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79600
Forward Sequence:	CTCTCGTCAACGCTGGACTT
Reverse Sequence:	TGGTGACTTCACCTGCATCTGTG
ACCN:	NM_001082538 , NM_001082538.1 , NM_001082538.2 , BC040113 , BC040113.1 , BC062611 , BC062611.1 , BC030993 , BC033811 , BC044885 , BP275624 , NM_001082538.3
UniProt ID:	Q2MV58
Synonyms:	JBTS13; TECT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)