

## Product datasheet for **HP230319**

### CCDC46 (CEP112) Human qPCR Primer Pair (NM\_145036)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	201134
Forward Sequence:	ATGGCAACTGGAGGACATCAGAC
Reverse Sequence:	GGCTATCCTTCTCTGCGTTATGC
ACCN:	<a href="#">BC029524</a> , <a href="#">NM_145036</a> , <a href="#">NM_145036.1</a> , <a href="#">NM_145036.2</a> , <a href="#">BC029524.1</a> , <a href="#">BC037991</a> , <a href="#">BC131805</a> , <a href="#">BG203878</a> , <a href="#">BG531260</a> , <a href="#">BX097629</a>
UniProt ID:	<a href="#">Q8N8E3</a>
Synonyms:	CCDC46; MACOCO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)