

Product datasheet for **HP230359**

ABR Human qPCR Primer Pair (NM_021962)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29
Forward Sequence:	TCTCCAGGACTCTATGGCTTCC
Reverse Sequence:	CACCCTGGTTTTGGCTTTGCTG
ACCN:	NM_021962 , NM_021962.1 , NM_021962.2 , NM_021962.3 , NM_021962.4 , BC036807 , BC141909 , BC156483 , BC172517 , BE675741 , BM547076 , BM827672 , NM_021962.5
UniProt ID:	Q12979
Synonyms:	MDB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)