

Product datasheet for **HP230375**

ALG3 Human qPCR Primer Pair (NM_005787)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10195
Forward Sequence:	CTTTGCTGTGCTCTACCTGGCT
Reverse Sequence:	CGCAGCACAAAGATGGAGTGGA
ACCN:	NM_005787 , NM_005787.1 , NM_005787.2 , NM_005787.3 , NM_005787.4 , NM_005787.5 , BC002839 , BC004313 , BC021062 , BM920198 , BX339649 , NM_005787.6
UniProt ID:	Q92685
Synonyms:	CDG1D; CDGS4; CDGS6; D16Ertd36e; not; Not56; NOT56L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)