

## Product datasheet for **HP230442**

### FAM48A (SUPT20H) Human qPCR Primer Pair (NM\_017569)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55578
Forward Sequence:	GACAGCCAGATGTCTCCATCAC
Reverse Sequence:	GAGCTGTGTGACATCTTGACCG
ACCN:	<u><a href="#">NM_017569</a></u> , <u><a href="#">NM_017569.1</a></u> , <u><a href="#">NM_017569.2</a></u> , <u><a href="#">NM_017569.3</a></u> , <u><a href="#">BC030686</a></u> , <u><a href="#">BC030686.2</a></u> , <u><a href="#">BC001145</a></u>
UniProt ID:	<u><a href="#">Q8NEM7</a></u>
Synonyms:	C13; C13orf19; FAM48A; FP757; P38IP; SPT20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)