

Product datasheet for **HP230470**

CDKL5 Human qPCR Primer Pair (NM_003159)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6792
Forward Sequence:	ATCCTGTGCCAAGGAGAGTTGG
Reverse Sequence:	GGTTGGTCCAGAACTGCTCTC
ACCN:	<u>NM_003159</u> , <u>NM_003159.1</u> , <u>NM_003159.2</u> , <u>BC010966</u> , <u>BC036091</u> , <u>BC152846</u>
UniProt ID:	<u>O76039</u>
Synonyms:	CFAP247; DEE2; EIEE2; ISSX; STK9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)