

Product datasheet for **HP230480**

CARD9 Human qPCR Primer Pair (NM_052813)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64170
Forward Sequence:	TCCGACCTGGAAGATGGCTCAC
Reverse Sequence:	CAGAGCTGCAAAGGGCTGTTTC
ACCN:	<u>NM_052813</u> , <u>NM_052813.1</u> , <u>NM_052813.2</u> , <u>NM_052813.3</u> , <u>NM_052813.4</u> , <u>BC008877</u> , <u>BC035805</u> , <u>BI524042</u> , <u>BM148244</u> , <u>NM_052813.5</u>
UniProt ID:	<u>Q9H257</u>
Synonyms:	CANDF2; hCARD9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)