

Product datasheet for **HP230498**

Glutaredoxin 1 (GLRX) Human qPCR Primer Pair (NM_002064)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2745
Forward Sequence:	ACAGCCACCAACCACACTAACG
Reverse Sequence:	GAGACTAGATCACTGCATCCGC
ACCN:	<u>BC010965</u> , <u>NM_002064</u> , <u>NM_002064.1</u> , <u>NM_002064.2</u> , <u>BC010965.1</u> , <u>BC005304</u> , <u>BC106075</u> , <u>BE738446</u> , <u>BG547098</u> , <u>BP196332</u> , <u>BT006689</u> , <u>BU949348</u> , <u>BX106369</u> , <u>NM_002064.3</u>
UniProt ID:	<u>P35754</u>
Synonyms:	GRX; GRX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)