

Product datasheet for **HP230517**

HIPK3 Human qPCR Primer Pair (NM_005734)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10114
Forward Sequence:	CTACAGATCCGACCAGGAGTTC
Reverse Sequence:	TGTGAACCAGCCCACTCTCAG
ACCN:	NM_005734 , NM_005734.1 , NM_005734.2 , NM_005734.3 , NM_005734.4 , BC030112 , BC146311 , BC156873 , BF332205 , BX117467 , NM_005734.5
UniProt ID:	Q9H422
Synonyms:	DYRK6; FIST3; PKY; YAK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)