

Product datasheet for **HP230527**

Kallikrein 6 (KLK6) Human qPCR Primer Pair (NM_002774)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5653
Forward Sequence:	GGTCCTTATCCATCCACTGTGG
Reverse Sequence:	GAACTCTCCCTTTGCCGAAGGT
ACCN:	<u>NM_002774</u> , <u>NM_002774.1</u> , <u>NM_002774.2</u> , <u>NM_002774.3</u> , <u>BC015525</u> , <u>BC015525.2</u> , <u>BM763657</u> , <u>BP253007</u> , <u>BT006852</u> , <u>NM_002774.4</u>
UniProt ID:	<u>Q92876</u>
Synonyms:	Bssp; hK6; Klk7; PRSS9; PRSS18; SP59
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)