

Product datasheet for **HP230529**

OCIAD1 Human qPCR Primer Pair (NM_017830)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54940
Forward Sequence:	CTGTGCCTTTGGCTGCAACAAG
Reverse Sequence:	AGTTTTCCAGCAAAGTATCCCATG
ACCN:	NM_017830 , NM_017830.1 , NM_017830.2 , NM_017830.3 , BC003409 , BC003409.1 , BC088361 , BI549536 , BM819316 , BM999322 , BQ272551 , NM_017830.4
UniProt ID:	Q9NX40
Synonyms:	ASRIJ; OCIA; TPA018
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)