

Product datasheet for **HP230534**

CLEC10A Human qPCR Primer Pair (NM_182906)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10462
Forward Sequence:	GGTGGATGGAACAGACTATGCG
Reverse Sequence:	GGATGGAAGTGAGCACAGTCCT
ACCN:	<u>NM_182906</u> , <u>NM_182906.1</u> , <u>NM_182906.2</u> , <u>NM_182906.3</u> , <u>BC039011</u> , <u>BC027858</u> , <u>BI910082</u> , <u>BM974488</u> , <u>BU687211</u>
UniProt ID:	<u>Q8IUN9</u>
Synonyms:	CD301; CLECSF13; CLECSF14; HML; HML2; MGL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)