

Product datasheet for **HP230567**

DIS3L Human qPCR Primer Pair (NM_001143688)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115752
Forward Sequence:	AGCAGAAACCTCCAGGACTTC
Reverse Sequence:	CTGTTTTCCACCAGGATGGTTGC
ACCN:	<u>NM_001143688</u> , <u>NM_001143688.1</u> , <u>NM_001143688.2</u> , <u>BC022089</u> , <u>BC022089.2</u> , <u>BC014124</u> , <u>BG722405</u> , <u>BM833083</u> , <u>BP219623</u> , <u>NM_001143688.3</u>
UniProt ID:	<u>Q8TF46</u>
Synonyms:	DIS3L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)