

Product datasheet for **HP230568**

FBXO11 Human qPCR Primer Pair (NM_025133)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80204
Forward Sequence:	ATCATGGACGTGATGTTGGTG
Reverse Sequence:	CCACTGTAGGGTTAGCATAGGC
ACCN:	<u>NM_025133</u> , <u>NM_025133.1</u> , <u>NM_025133.2</u> , <u>NM_025133.3</u> , <u>NM_025133.4</u> , <u>BC012728</u> , <u>BC043258</u> , <u>BC130445</u> , <u>BC136480</u>
UniProt ID:	<u>Q86XK2</u>
Synonyms:	FBX11; IDDFBA; PRMT9; UBR6; UG063H01; VIT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)