

Product datasheet for **HP230577**

Neuroserpin (SERPINI1) Human qPCR Primer Pair (NM_005025)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5274
Forward Sequence:	CTACCCAGAAAGAAATCCGCCAC
Reverse Sequence:	TGGCTCTCTTTAGCAGTTACCATG
ACCN:	NM_005025 , NM_005025.1 , NM_005025.2 , NM_005025.3 , NM_005025.4 , BC018043 , BC018043.1 , BQ086144
UniProt ID:	Q99574
Synonyms:	HNS-S1; HNS-S2; neuroserpin; PI12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)