

Product datasheet for **HP230751**

HIP55 (DBNL) Human qPCR Primer Pair (NM_014063)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28988
Forward Sequence:	AGGATGTGGAGCCTGAGTGCAT
Reverse Sequence:	GACACGGCATTGGTCTTCTGGT
ACCN:	NM_014063 , NM_014063.1 , NM_014063.2 , NM_014063.3 , NM_014063.4 , NM_014063.5 , NM_014063.6 , BC011677 , BC023589 , BC031687 , NM_014063.7
UniProt ID:	Q9UJU6
Synonyms:	ABP1; HIP-55; HIP55; SH3P7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)