

Product datasheet for **HP230767**

CD8B Human qPCR Primer Pair (NM_172099)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	926
Forward Sequence:	AAGACAGTGGCATCTACTTCTGC
Reverse Sequence:	TTGAGGGTGGACTTCTTGGTGG
ACCN:	<u>NM_172099</u> , <u>NM_172099.1</u> , <u>NM_172099.2</u> , <u>BC100911</u> , <u>BC100912</u> , <u>BC100913</u> , <u>BC100914</u> , <u>BP392114</u> , <u>BQ004743</u>
UniProt ID:	<u>P10966</u>
Synonyms:	CD8B1; Leu2; Ly3; LYT3; MGC119115
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)