

Product datasheet for **HP230769**

PARP3 Human qPCR Primer Pair (NM_005485)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10039
Forward Sequence:	TCTCTGAGCAGGAGAAGACGGT
Reverse Sequence:	TGTGGTTGCTGCCAGTCTGTTC
ACCN:	NM_005485 , NM_005485.1 , NM_005485.2 , NM_005485.3 , NM_005485.4 , NM_005485.5 , BC014260 , BF375394 , BG818516 , BI554046 , BU632340
UniProt ID:	Q9Y6F1
Synonyms:	ADPRT3; ADPRTL2; ADPRTL3; ARTD3; IRT1; PADPRT-3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)