

Product datasheet for **HP230776**

PRMT5 Human qPCR Primer Pair (NM_006109)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10419
Forward Sequence:	CTAGACCGAGTACCAGAAGAGG
Reverse Sequence:	CAGCATACAGCTTTATCCGCCG
ACCN:	NM_006109 , NM_006109.1 , NM_006109.2 , NM_006109.3 , NM_006109.4 , BC025979 , BC025979.2 , BC005820 , NM_006109.5
UniProt ID:	O14744
Synonyms:	HRMT1L5; HSL7; IBP72; JBP1; SKB1; SKB1Hs
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)