

Product datasheet for **HP230781**

Glucosidase 2 subunit beta (PRKCSH) Human qPCR Primer Pair (NM_002743)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5589
Forward Sequence:	GCGTCATCTGTGAGAACACCTG
Reverse Sequence:	CGTGCCTTCTCCAGTCCTCAA
ACCN:	<u>NM_002743</u> , <u>NM_002743.1</u> , <u>NM_002743.2</u> , <u>BC013586</u> , <u>BC015154</u> , <u>BC128412</u> , <u>BF984334</u> , <u>BG323695</u> , <u>BI254784</u> , <u>BP242217</u> , <u>BT009858</u>
UniProt ID:	<u>P14314</u>
Synonyms:	AGE-R2; G19P1; GIIB; PCLD; PCLD1; PKCSH; PLD1; VASAP-60
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)