

Product datasheet for **HP230782**

Zyxin (ZYG) Human qPCR Primer Pair (NM_003461)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7791
Forward Sequence:	TTCCACATCGCCTGCTTACCT
Reverse Sequence:	CGCAGGTGTACACTTCTCCAG
ACCN:	<u>NM_003461</u> , <u>NM_003461.1</u> , <u>NM_003461.2</u> , <u>NM_003461.3</u> , <u>NM_003461.4</u> , <u>BC008743</u> , <u>BC002323</u> , <u>BC009360</u> , <u>BC010031</u> , <u>BC017183</u> , <u>BU844538</u> , <u>NM_003461.5</u>
UniProt ID:	<u>Q15942</u>
Synonyms:	ESP-2; HED-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)