

Product datasheet for **HP230798**

GAB2 Human qPCR Primer Pair (NM_080491)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9846
Forward Sequence:	TCAGCAGAGACCGCCAATCAGT
Reverse Sequence:	GGTACTCGTAGGTCTCACAGGA
ACCN:	<u>NM_080491</u> , <u>NM_080491.1</u> , <u>NM_080491.2</u> , <u>BC131711</u> , <u>BC152459</u> , <u>BY799914</u> , <u>NM_080491.3</u>
UniProt ID:	<u>Q9UQC2</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)