

Product datasheet for **HP230828**

U2AF1L3 (U2AF1L4) Human qPCR Primer Pair (NM_144987)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	199746
Forward Sequence:	GTGGCTGAACTCAGTAACCGCT
Reverse Sequence:	CATAGAGCTGCCTCTGGAGGTT
ACCN:	<u>NM_144987</u> , <u>NM_144987.1</u> , <u>NM_144987.2</u> , <u>NM_144987.3</u> , <u>BC021186</u> , <u>BC021186.1</u> , <u>BC010865</u> , <u>BC049208</u> , <u>BC062474</u> , <u>BM696851</u>
UniProt ID:	<u>Q8WU68</u>
Synonyms:	U2AF1-RS3; U2AF1L3; U2AF1L3V1; U2AF1RS3; U2af26
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)