

Product datasheet for **HP230871**

ZNF331 Human qPCR Primer Pair (NM_018555)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55422
Forward Sequence:	ACTCTGCTCAGAGGGACCTGTA
Reverse Sequence:	TCCTTTTGGAAGCCCTTATTTTCATG
ACCN:	NM_018555 , NM_018555.1 , NM_018555.2 , NM_018555.3 , NM_018555.4 , NM_018555.5 , BC009433 , BC009433.1 , BC009365 , BF057584 , BP369529 , BQ003434 , NM_018555.6
UniProt ID:	Q9NQX6
Synonyms:	RITA; ZNF361; ZNF463
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)