

## Product datasheet for **HP230879**

### Neurofilament (NEFM) Human qPCR Primer Pair (NM\_005382)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4741
Forward Sequence:	ACAACCACGACCTCAGCAGCTA
Reverse Sequence:	GTTGAGGAGGTCCTGGTATTCG
ACCN:	<a href="#">NM_005382</a> , <a href="#">NM_005382.1</a> , <a href="#">NM_005382.2</a> , <a href="#">BC096757</a> , <a href="#">BC002421</a> , <a href="#">BC071752</a> , <a href="#">BI496345</a> , <a href="#">BP247321</a> , <a href="#">BP311096</a>
UniProt ID:	<a href="#">P07197</a>
Synonyms:	NEF3; NF-M; NFM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)