

## Product datasheet for **HP230883**

### NCBP2 Human qPCR Primer Pair (NM\_007362)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22916
Forward Sequence:	CCGAATCATTGCGACAGACTGG
Reverse Sequence:	CTGTGCCAGTTTTCCATAGCCTC
ACCN:	<a href="#">NM_007362</a> , <a href="#">NM_007362.1</a> , <a href="#">NM_007362.2</a> , <a href="#">NM_007362.3</a> , <a href="#">BC001255</a> , <a href="#">BC001255.1</a> , <a href="#">BP248419</a> , <a href="#">BT006842</a> , <a href="#">BX537372</a> , <a href="#">BX647994</a> , <a href="#">NM_007362.5</a>
UniProt ID:	<a href="#">P52298</a>
Synonyms:	CBC2; CBP20; NIP1; PIG55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)