

Product datasheet for **HP230896**

PCK2 Human qPCR Primer Pair (NM_004563)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5106
Forward Sequence:	TAGTGCCTGTGGCAAGACCAAC
Reverse Sequence:	GAAGCCGTTCTCAGGGTTGATG
ACCN:	<u>NM_004563</u> , <u>NM_004563.1</u> , <u>NM_004563.2</u> , <u>NM_004563.3</u> , <u>BC001454</u> , <u>BQ957219</u> , <u>BU153045</u> , <u>BU623248</u> , <u>BX248272</u> , <u>NM_004563.4</u>
UniProt ID:	<u>Q16822</u>
Synonyms:	PEPCK; PEPCK-M; PEPCK2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)