

## Product datasheet for **HP230905**

### Prostatic Acid Phosphatase (ACPP) Human qPCR Primer Pair (NM\_001134194)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55
Forward Sequence:	CAAGACTGGTCCACGGAGTGTA
Reverse Sequence:	AGCAGAGTCCACGGCGAATGTG
ACCN:	<a href="#">NM_001134194</a> , <a href="#">NM_001134194.1</a> , <a href="#">BC016344</a> , <a href="#">BC016344.1</a> , <a href="#">BC007460</a> , <a href="#">BC008493</a> , <a href="#">BC017877</a> , <a href="#">BQ100881</a> , <a href="#">NM_001134194.2</a>
UniProt ID:	<a href="#">P15309</a>
Synonyms:	5'-NT; ACP-3; ACP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)