

## Product datasheet for **HP230945**

### **L3MBTL3 Human qPCR Primer Pair (NM\_032438)**

#### **Product data:**

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 84456   |
| Forward Sequence: | TTCGCAGAGAGCACGGAGGAAA  |
| Reverse Sequence: | ACCGCTTTCTCCTCTCCAGGT   |
| ACCN:             | <u>NM_032438</u> , <u>NM_032438.1</u> , <u>NM_032438.2</u> , <u>NM_032438.3</u> , <u>BC060845</u> , <u>BC060845.1</u> ,<br><u>NM_032438.4</u>   |
| UniProt ID:       | <u>Q96JM7</u>   |
| Synonyms:         | MBT-1; MBT1   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)