

Product datasheet for **HP230983**

PTP4A3 Human qPCR Primer Pair (NM_032611)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11156
Forward Sequence:	GACCTATGACAAAACGCCGCTG
Reverse Sequence:	CAGCCAGTCTTCCACTACCTTG
ACCN:	BC003105 , NM_032611 , NM_032611.1 , BC003105.1 , BE778111 , BQ637038 , BT007303
UniProt ID:	O75365
Synonyms:	PRL-3; PRL-R; PRL3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)