

Product datasheet for **HP231016**

p38 (CRK) Human qPCR Primer Pair (NM_016823)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1398
Forward Sequence:	CCAATGCCTACGACAAGACAGC
Reverse Sequence:	TGGGAAGTGACCTCGTTTGCCA
ACCN:	NM_016823 , NM_016823.1 , NM_016823.2 , NM_016823.3 , BC008506 , BC008506.1 , BC009837 , BC009837.1 , BC001718 , BT007277 , NM_016823.4
UniProt ID:	P46108
Synonyms:	CRKII; p38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)