

Product datasheet for **HP231021**

MGAT2 Human qPCR Primer Pair (NM_002408)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4247
Forward Sequence:	CCAGGGAATTGACAACGTCCTC
Reverse Sequence:	AGGGTCACTACCTGGAAACTCG
ACCN:	<u>NM_002408</u> , <u>NM_002408.1</u> , <u>NM_002408.2</u> , <u>NM_002408.3</u> , <u>BC006390</u> , <u>BC006390.2</u> , <u>BC018912</u> , <u>BI764270</u> , <u>BQ023756</u> , <u>NM_002408.4</u>
UniProt ID:	<u>Q10469</u>
Synonyms:	CDG2A; CDGS2; GLCNACTII; GNT-II; GNT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)