

## Product datasheet for **HP231057**

### Glutathione Transferase zeta 1 (GSTZ1) Human qPCR Primer Pair (NM\_145870)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2954
Forward Sequence:	AGGACTTCCAGGCACTGAATCC
Reverse Sequence:	CTTTGGTCCTGAGGCAGAAGT
ACCN:	<a href="#">NM_145870</a> , <a href="#">NM_145870.1</a> , <a href="#">NM_145870.2</a> , <a href="#">BC001453</a> , <a href="#">BC001453.1</a> , <a href="#">BC052245</a> , <a href="#">BC069047</a> , <a href="#">BE958109</a> , <a href="#">BG330365</a> , <a href="#">BG772450</a> , <a href="#">BI559916</a> , <a href="#">BI829992</a> , <a href="#">NM_145870.3</a>
UniProt ID:	<a href="#">O43708</a>
Synonyms:	GSTZ1-1; MAAI; MAAID; MAI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)