

## Product datasheet for **HP231065**

### Centlein (CNTLN) Human qPCR Primer Pair (NM\_017738)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54875
Forward Sequence:	GGTCGTCACCACACTGTTCTCA
Reverse Sequence:	CCAAAGTTTGAGAGGTCTGTGCT
ACCN:	<a href="#">NM_017738</a> , <a href="#">NM_017738.1</a> , <a href="#">NM_017738.2</a> , <a href="#">NM_017738.3</a> , <a href="#">BC156251</a> , <a href="#">BG575335</a> , <a href="#">BU075047</a> , <a href="#">BU162362</a> , <a href="#">BU663791</a> , <a href="#">BX647069</a> , <a href="#">NM_017738.4</a>
UniProt ID:	<a href="#">Q9NXG0</a>
Synonyms:	bA340N12.1; C9orf39; C9orf101
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)