

Product datasheet for **HP231082**

Liver Carboxylesterase 1 (CES1) Human qPCR Primer Pair (NM_001266)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1066
Forward Sequence:	ATCCAAGGCAACT
Reverse Sequence:	GACAGTGTCTGTTCTCCT
ACCN:	NM_001266 , NM_001266.1 , NM_001266.2 , NM_001266.3 , NM_001266.4 , BC012418 , BC012418.1 , BC009706 , BC110338 , NM_001266.5
UniProt ID:	P23141
Synonyms:	ACAT; CE-1; CEH; CES2; hCE-1; HMSE; HMSE1; PCE-1; REH; SES1; TGH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)