

Product datasheet for **HP231125**

Plakophilin 2 (PKP2) Human qPCR Primer Pair (NM_004572)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5318
Forward Sequence:	CAGGTGCTGAAGCAAACCAGAG
Reverse Sequence:	CCAAAGTGCTGGGATTACAGGC
ACCN:	NM_004572 , NM_004572.1 , NM_004572.2 , NM_004572.3 , BC094762 , BC094762.1 , BC070083 , BC126199 , BC143966 , BC143967 , BX956935 , NM_004572.4
UniProt ID:	Q99959
Synonyms:	ARVD9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)