

Product datasheet for **HP231168**

Fibulin 1 (FBLN1) Human qPCR Primer Pair (NM_006486)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2192
Forward Sequence:	GTGGCTACCATCTCAACGAGGA
Reverse Sequence:	CTTGCAATTCGCAGCGAAACTG
ACCN:	NM_006486 , NM_006486.1 , NM_006486.2 , BC022497 , BC022497.1 , BX403501 , BX424750 , BX443677 , NM_006486.3
UniProt ID:	P23142
Synonyms:	FBLN; FIBL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)