

Product datasheet for **HP231267**

SSBP3 Human qPCR Primer Pair (NM_018070)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23648
Forward Sequence:	TCAGGCGACATAGACGGACTTC
Reverse Sequence:	CACACTCATCGTCATGCTTGGAG
ACCN:	NM_018070 , NM_018070.2 , NM_018070.3 , NM_018070.4 , BC003605 , BC003605.1 , BC066365 , BC066365.1 , BC004335 , BC127866 , BF060746 , BQ051652 , NM_018070.5
UniProt ID:	Q9BWW4
Synonyms:	CSDP; SSDP; SSDP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)