

## Product datasheet for **HP231346**

### **SRC3 (NCOA3) Human qPCR Primer Pair (NM\_181659)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	8202
<b>Forward Sequence:</b>	GGACTAAGCAACAGGTGTTTCAAG
<b>Reverse Sequence:</b>	ACTGGAGGACTTGAGCCAACAG
<b>ACCN:</b>	<u><a href="#">NM_181659</a></u> , <u><a href="#">NM_181659.1</a></u> , <u><a href="#">NM_181659.2</a></u> , <u><a href="#">BC092516</a></u> , <u><a href="#">BC119001</a></u> , <u><a href="#">BC122547</a></u> , <u><a href="#">BM692321</a></u> , <u><a href="#">BU588121</a></u> , <u><a href="#">NM_181659.3</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q9Y6Q9</a></u>
<b>Synonyms:</b>	ACTR; AIB-1; AIB1; bHLHe42; CAGH16; CTG26; KAT13B; pCIP; RAC3; SRC-3; SRC3; TNRC14; TNRC16; TRAM-1
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)