

Product datasheet for **HP231371**

RAB5IF Human qPCR Primer Pair (NM_199483)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55969
Forward Sequence:	AGTCTCCGTCTGGAGTAAGGTG
Reverse Sequence:	CCTGCTATTCCCAAGAACCCTC
ACCN:	<u>NM_199483</u> , <u>NM_199483.1</u> , <u>NM_199483.2</u> , <u>BC001871</u> , <u>BC004446</u> , <u>BG395753</u> , <u>BI223439</u> , <u>BQ086248</u> , <u>BU174619</u> , <u>BU189639</u> , <u>BU194101</u> , <u>BU848918</u> , <u>NM_199483.3</u>
UniProt ID:	<u>Q9BUV8</u>
Synonyms:	C20orf24; PNAS-11; RCAF1; RIP5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)