

## Product datasheet for **HP231378**

### STAT3 Human qPCR Primer Pair (NM\_139276)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6774
Forward Sequence:	CTTTGAGACCGAGGTGTATCACC
Reverse Sequence:	GGTCAGCATGTTGTACCACAGG
ACCN:	<a href="#">NM_139276</a> , <a href="#">NM_139276.1</a> , <a href="#">NM_139276.2</a> , <a href="#">BC014482</a> , <a href="#">BC014482.1</a> , <a href="#">BC000627</a> , <a href="#">BC008044</a> , <a href="#">BC029783</a> , <a href="#">BC067119</a> , <a href="#">BC107775</a> , <a href="#">BI461226</a> , <a href="#">NM_139276.3</a>
UniProt ID:	<a href="#">P40763</a>
Synonyms:	ADMIO; ADMIO1; APRF; HIES
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)