

Product datasheet for **HP231384**

OGG1 Human qPCR Primer Pair (NM_016821)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4968
Forward Sequence:	GGCTCAACTGTATCACCCTGG
Reverse Sequence:	GGCGATGTTGTTGTTGGAGGAAC
ACCN:	<u>NM_016821</u> , <u>NM_016821.1</u> , <u>NM_016821.2</u> , <u>BC000657</u> , <u>BC000657.2</u> , <u>BM847736</u>
UniProt ID:	<u>O15527</u>
Synonyms:	HMMH; HOGG1; MUTM; OGH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)