

Product datasheet for **HP231402**

FNTA Human qPCR Primer Pair (NM_002027)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2339
Forward Sequence:	AGCAGCACAAAGGAAGAGATGGC
Reverse Sequence:	ATCAGCCCATTCTGCTCTGTCC
ACCN:	NM_002027 , NM_002027.1 , NM_002027.2 , BC084566 , BC084566.1 , BC017029 , BC037295 , BI913631 , BM686906 , BT009854 , BX648381 , NM_002027.3
UniProt ID:	P49354
Synonyms:	FPTA; PGGT1A; PTAR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)