

Product datasheet for **HP231406**

SH3KBP1 Human qPCR Primer Pair (NM_031892)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	30011
Forward Sequence:	GCAGTTCGCTATCTGGCATCCT
Reverse Sequence:	GTCTGCTTGTTGGTCGGATGACT
ACCN:	NM_031892 , NM_031892.1 , NM_031892.2 , BC015806 , BC015806.1 , BC041055 , BC050663 , BP224209 , NM_031892.3
UniProt ID:	Q96B97
Synonyms:	AGMX2; CD2BP3; CIN85; GIG10; HSB-1; HSB1; IMD61; MIG18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)