

## Product datasheet for **HP231409**

### hnRNP A2B1 (HNRNPA2B1) Human qPCR Primer Pair (NM\_031243)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3181
Forward Sequence:	CAGCAACCTTCTAACTACGGTCC
Reverse Sequence:	CACTGCCTCCTGGACCATAGTT
ACCN:	<u><a href="#">NM_031243</a></u> , <u><a href="#">NM_031243.1</a></u> , <u><a href="#">NM_031243.2</a></u> , <u><a href="#">BC000506</a></u> , <u><a href="#">BC045724</a></u> , <u><a href="#">BX537494</a></u> , <u><a href="#">BX538099</a></u> , <u><a href="#">BX641084</a></u> , <u><a href="#">NM_031243.3</a></u>
UniProt ID:	<u><a href="#">P22626</a></u>
Synonyms:	HNRNPA2; HNRNPB1; HNRPA2; HNRPA2B1; HNRPB1; IBMPFD2; RNPA2; SNRPB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)