

Product datasheet for **HP231424**

MACROD2 Human qPCR Primer Pair (NM_080676)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140733
Forward Sequence:	GAAGATGCCGTGATTGTGGAGC
Reverse Sequence:	TGGAGCCTTCACTTTTCCAGG
ACCN:	BC018687 , NM_080676 , NM_080676.2 , NM_080676.3 , NM_080676.4 , NM_080676.5 , BC018687.1 , BC035876 , BC036403 , BC101216 , BC101217 , BC101218 , BC101219 , BC126936 , BC128036 , BX281637
UniProt ID:	A1Z1Q3
Synonyms:	C2orf133; C20orf133
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)