

## Product datasheet for **HP231456**

### Vinculin (VCL) Human qPCR Primer Pair (NM\_014000)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7414
Forward Sequence:	TGAGCAAGCACAGCGGTGGATT
Reverse Sequence:	TCGGTCACACTTGCGGAGAAGA
ACCN:	<u><a href="#">NM_014000</a></u> , <u><a href="#">NM_014000.1</a></u> , <u><a href="#">NM_014000.2</a></u> , <u><a href="#">BC039174</a></u> , <u><a href="#">BC039174.1</a></u> , <u><a href="#">BC013847</a></u> , <u><a href="#">BX537994</a></u> , <u><a href="#">NM_014000.3</a></u>
UniProt ID:	<u><a href="#">P18206</a></u>
Synonyms:	CMD1W; CMH15; HEL114; MV; MVCL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)