

## Product datasheet for **HP231519**

### TCP1 alpha (TCP1) Human qPCR Primer Pair (NM\_030752)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6950
Forward Sequence:	TTTGTGGAGGCTGGTGCTATGG
Reverse Sequence:	CTTCTGCCTGTCCCAACATTGC
ACCN:	<a href="#">BC000665</a> , <a href="#">NM_030752</a> , <a href="#">NM_030752.1</a> , <a href="#">NM_030752.2</a> , <a href="#">BC000665.1</a> , <a href="#">BM697304</a> , <a href="#">BM827784</a> , <a href="#">BT006969</a> , <a href="#">NM_030752.3</a>
UniProt ID:	<a href="#">P17987</a>
Synonyms:	CCT-alpha; CCT1; CCTa; D6S230E; TCP-1-alpha
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)