

Product datasheet for **HP231539**

APPBP1 (NAE1) Human qPCR Primer Pair (NM_003905)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8883
Forward Sequence:	CCTGTTTCGAGGCACAATTCCTG
Reverse Sequence:	TTACCCACAGCGGCAGCATCTT
ACCN:	NM_003905 , NM_003905.1 , NM_003905.2 , NM_003905.3 , BC013301 , BC013301.2 , BC000480 , BC041323 , BM045535 , BM785899 , NM_003905.4
UniProt ID:	Q13564
Synonyms:	A-116A10.1; APPBP1; HPP1; ula-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)