

Product datasheet for **HP231644**

MCK10 (DDR1) Human qPCR Primer Pair (NM_013994)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	780
Forward Sequence:	GATCTCGACTCCGCTTCAAGGA
Reverse Sequence:	CAAAGGGTGCCCTTACGCACA
ACCN:	<u>NM_013994</u> , <u>NM_013994.1</u> , <u>NM_013994.2</u> , <u>BC013400</u> , <u>BC013400.2</u> , <u>BC008716</u> , <u>BC070070</u> , <u>BI193181</u> , <u>BM875592</u>
UniProt ID:	<u>Q08345</u>
Synonyms:	CAK; CD167; DDR; EDDR1; HGK2; MCK10; NEP; NTRK4; PTK3; PTK3A; RTK6; TRKE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)